

Suppl. Figure 3 GC-MS data for dihydronepetalactone (4), 2 July 2014.

Analyzed on an HP 6890N GC coupled in series with an HP 5973 mass selective detector using a 30m DB-5 capillary column (250 μm x 0.25 μm film Thickness; Agilent Technologies, Wilmington, DE, USA), 50 $^{\circ}\text{C}$ for 5 min, to 280 $^{\circ}\text{C}$ at 10 $^{\circ}\text{C}/\text{min}$, hold 3 min.

